
3/97-Sf-

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This apparatus illustrates the Debye-Scherrer process used in the electron beam diffraction tube (555 17 or 555 234) with the aid of visible light.

Using an iris diaphragm, a narrow beam is screened out from white light emitted from a light source and made approximately parallel using a lens system (or one lens). After using a color filter to obtain a monochromatic light, the beam passes through a rotating lattice grating. The rotation blurs the diffraction image produced by a lattice grating – light spots arranged in a network pattern around the undiffracted central beam – to form two concentric rings arranged around the central light spot.

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Durchscheinender Schirm 441 53

1 Description, scope of supply

- ① Lattice grating (dia. 3 mm) in ball bearing-mounted flywheel (dia. 10 cm)
- ② Iris diaphragm (dia. 1 mm) in frame (5 cm x 5 cm)
- ③ 2 color filters (8 cm x 10 cm), plastic

2 Operation

Additionally required:

Holder for lattice grating ①
either (Fig. 2.1)

Stand for electron tubes 555 05

or (Fig. 2.2)

Holder with spring clips 460 22

High-intensity illumination apparatus with imaging system

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Carry out the experiment in a darkened room.
 Before placing the lattice grating ① in the beam path, focus the
 iris diaphragm ② sharply on the screen.

First observe the diffraction image of the stationary lattice gra-
 ting, and then the image of the rotating lattice grating.

Use both white and filtered light. Observe how the filter color
 influences the radius of the diffraction rings.

Analogies:

	Electron beam diffraction (with 555 17 or 234)	Optical analogon (with 555 19)
Beam	Hot electrons	Visible light, white
Focusing	Electron-optical	Lenses
Diaphragm	Iris diaphragm	Iris diaphragm
Monochromaticity	Wavelength λ determined by anode voltage U : $\lambda = f(U)$	color filter
Diffraction grating	Polycrystalline graphite	Rotating lattice grating
Diffraction pattern	Concentric circles for which the diameter changes acc. to $\lambda = f(U)$	Concentric circles for which the diameter changes with the color filter


Fig. 2.1


Fig. 2.2